

Product Datasheet

Mouse Enkephalin (ENK) ELISA Kit (orb1666600)

Description	Mouse Enkephalin (ENK) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 93.7-6000 pg/mL. Sensitivity: 25.35 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.