

Product Datasheet

Human Esophageal Cancer Related Gene 4 (ECRG4) ELISA Kit (orb1666668)

Description	Human Esophageal Cancer Related Gene 4 (ECRG4) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 30 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.