

Product Datasheet

Guinea pig Eosinophil Chemotactic Factor (ECF) ELISA Kit (orb1666673)

Description	Guinea pig Eosinophil Chemotactic Factor (ECF) ELISA Kit
Reactivity	Guinea pig
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.2 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.