

Product Datasheet

CMKLR1 antibody (orb166668)

Description	Rabbit polyclonal antibody to CMKLR1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CMKLR1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42322
Uniprot ID	Q99788
NCBI	NP_001135816.1 , NP_001135815.1 , NP_001135817.1 , NP_004063.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

