

Product Datasheet

Human Dolichyl Phosphate-N-Acetylglucosaminephosphotransferase 1 (DPAGT1) ELISA Kit (orb1666724)

Description Human Dolichyl Phosphate-N-Acetylglucosaminephosphotransferase 1 (DPAGT1) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL.
Sensitivity: 0.059 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.

Expiration Date Please enquire.