

Product Datasheet

RNASE1 antibody (orb166675)

Description	Rabbit polyclonal antibody to RNASE1 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RNASE1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	17644
Uniprot ID	P07998
NCBI	NP_937877.1 , NP_937878.1 , NP_002924.1 , NP_937875.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

