

Product Datasheet

Human Defensin Beta 103A (DEFb103A) ELISA Kit (orb1666801)

Description	Human Defensin Beta 103A (DEFb103A) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 46.8-3000 pg/mL. Sensitivity: 12.55 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.