

Product Datasheet

LRCH1 antibody (orb166689)

Description	Rabbit polyclonal antibody to LRCH1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LRCH1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	80875
Uniprot ID	Q9Y2L9
NCBI	NP_001157685.1 , NP_001157683.1 , NP_055931.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

