

Product Datasheet

MEG11 antibody (orb166690)

Description	Rabbit polyclonal antibody to MEG11 (N-term).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MEGF11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	110844
Uniprot ID	A6BM72
NCBI	NP_115821.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

