

Product Datasheet

WBSCR22 antibody (orb166693)

Description	Rabbit polyclonal antibody to WBSCR22 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	BUD23
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31880
Uniprot ID	O43709
NCBI	NP_001189489.1 , NP_059998.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

