

Product Datasheet

Bovine Casein Kappa (CSN3) ELISA Kit (orb1666972)

Description	Bovine Casein Kappa (CSN3) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 µg/mL. Sensitivity: 0.101 µg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: milk.
Expiration Date	Please enquire.