

Product Datasheet

Mouse Cartilage Associated Protein (CRTAP) ELISA Kit (orb1666988)

Description	Mouse Cartilage Associated Protein (CRTAP) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.056 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.