

Product Datasheet

Human COP9 Constitutive Photomorphogenic Homolog Subunit 2 (COPS2) ELISA Kit (orb1667080)

Description	Human COP9 Constitutive Photomorphogenic Homolog Subunit 2 (COPS2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.065 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.