

Product Datasheet

Rat Collagen Type IX Alpha 1 (COL9a1) ELISA Kit (orb1667089)

Description	Rat Collagen Type IX Alpha 1 (COL9a1) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.52 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.