

Product Datasheet

Rat Collagen Type I Alpha 2 (COL1a2) ELISA Kit (orb1667107)

Description	Rat Collagen Type I Alpha 2 (COL1a2) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.55 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.