

Product Datasheet

Rat Collagen Type XIX (COL19) ELISA Kit (orb1667115)

Description	Rat Collagen Type XIX (COL19) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.116 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.