

Product Datasheet

EDC3 antibody (orb166712)

Description	Rabbit polyclonal antibody to EDC3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	EDC3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56078
Uniprot ID	Q96F86
NCBI	NP_001135916.1 , NP_079359.2 , NP_001135915.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

