

Product Datasheet

KCNMB2 antibody (orb166713)

Description	Rabbit polyclonal antibody to KCNMB2 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KCNMB2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	27130
Uniprot ID	Q9Y691
NCBI	NP_001265840.1 , NP_852006.1 , NP_005823.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

