

Product Datasheet

Chicken Ciliary Neurotrophic Factor (CNTF) ELISA Kit (orb1667135)

Description	Chicken Ciliary Neurotrophic Factor (CNTF) ELISA Kit
Reactivity	Gallus
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5-4000 pg/mL. Sensitivity: 29.4 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.