

Product Datasheet

Human Cholinergic Receptor, Nicotinic, Beta 2 (CHRNb2) ELISA Kit (orb1667200)

Description	Human Cholinergic Receptor, Nicotinic, Beta 2 (CHRNb2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.054 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.