

## Product Datasheet

### Rat Chromogranin B (CHGB) ELISA Kit (orb1667216)

|                            |                                                                                                                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rat Chromogranin B (CHGB) ELISA Kit                                                                                                                                                                                         |
| <b>Reactivity</b>          | Rat                                                                                                                                                                                                                         |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                       |
| <b>Note</b>                | For research use only                                                                                                                                                                                                       |
| <b>Application notes</b>   | Detect Range: 0.625-40 ng/mL.<br>Sensitivity: 0.266 ng/mL. Assay Time:<br>1-4.5 h. Detection wavelength: 450 nm.<br>Sample Volume: 100 µL. Sample Type:<br>serum, plasma, tissue homogenates or<br>other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                                             |