

Product Datasheet

Mouse Transforming growth factor beta-1 ELISA Kit (orb1667230)

Description	Mouse Transforming growth factor beta-1 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 T. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.