

Product Datasheet

Human Cholesterol-25-Hydroxylase (CH25H) ELISA Kit (orb1667232)

Description Human Cholesterol-25-Hydroxylase (CH25H) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.625-40 ng/mL. Sensitivity: 0.252 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.

Expiration Date Please enquire.