

Product Datasheet

Human Calcitonin Gene Related Peptide (CGRP) ELISA Kit (orb1667235)

Description Human Calcitonin Gene Related Peptide (CGRP) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 15.6-1000 pg/mL. Sensitivity: 4.26 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Expiration Date Please enquire.