

Product Datasheet

BAD (phospho-Y113) antibody (orb166729)

Description	Rabbit polyclonal antibody to BAD (phospho-Y113)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	Bad
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22080
Uniprot ID	Q61337
NCBI	NP_031548.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

