

Product Datasheet

Bovine Cluster Of Differentiation 40 Ligand (CD40L) ELISA Kit (orb1667331)

Description	Bovine Cluster Of Differentiation 40 Ligand (CD40L) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.18 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.