

Product Datasheet

Human CD200 Receptor 2 (CD200R2) ELISA Kit (orb1667343)

Description	Human CD200 Receptor 2 (CD200R2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.781-50 ng/mL. Sensitivity: 0.34 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.