

Product Datasheet

Goat Caveolin 1 (CAV1) ELISA Kit (orb1667425)

Description	Goat Caveolin 1 (CAV1) ELISA Kit
Reactivity	Goat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.06 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.