

Product Datasheet

ERBB2 (phospho-Y1223) antibody (orb166745)

Description	Rabbit polyclonal antibody to ERBB2 (phospho-Y1223).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	ErbB2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	138579
Uniprot ID	P70424
NCBI	NP_001003817.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

