

Product Datasheet

Mouse Calpain 10 (CAPN10) ELISA Kit (orb1667471)

Description Mouse Calpain 10 (CAPN10) ELISA Kit

Reactivity Mouse

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL.
Sensitivity: 0.103 ng/mL. Assay Time:
1-4.5 h. Detection wavelength: 450
nm. Sample Volume: 100 µL. Sample
Type: tissue homogenates or other
biological fluids.

Expiration Date Please enquire.