

Product Datasheet

Mouse Calcium Channel, Voltage Dependent, L-Type, Alpha 1C Subunit (CACNa1C) ELISA Kit (orb1667500)

Description	Mouse Calcium Channel, Voltage Dependent, L-Type, Alpha 1C Subunit (CACNa1C) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.115 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.