

Product Datasheet

c-Kit (phospho-S739) antibody (orb166755)

Description	Rabbit polyclonal antibody to c-Kit (phospho-S739)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	Kit
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	109343
Uniprot ID	P05532
NCBI	NP_001116205.1 , NP_066922.2
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

