

Product Datasheet

Human Bone Marrow Stromal Cell Antigen 2 (BST2) ELISA Kit (orb1667577)

Description	Human Bone Marrow Stromal Cell Antigen 2 (BST2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.068 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.