

Product Datasheet

Bovine Bone Morphogenetic Protein Receptor 1A (BMPR1A) ELISA Kit (orb1667603)

Description	Bovine Bone Morphogenetic Protein Receptor 1A (BMPR1A) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.53 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.