

Product Datasheet

Goat Bone Morphogenetic Protein 4 (BMP4) ELISA Kit (orb1667615)

Description	Goat Bone Morphogenetic Protein 4 (BMP4) ELISA Kit
Reactivity	Goat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.148 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.