

Product Datasheet

Chicken Bone Morphogenetic Protein 4 (BMP4) ELISA Kit (orb1667616)

Description	Chicken Bone Morphogenetic Protein 4 (BMP4) ELISA Kit
Reactivity	Gallus
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 12.6 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.