

Product Datasheet

Mouse Alpha-2-Glycoprotein 1, Zinc Binding (aZGP1) ELISA Kit (orb1667707)

Description	Mouse Alpha-2-Glycoprotein 1, Zinc Binding (aZGP1) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.25-80 ng/mL. Sensitivity: 0.48 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.