

Product Datasheet

Mouse ATPase, H⁺/K⁺ Exchanging Beta Polypeptide (ATP4b) ELISA Kit (orb1667718)

Description	Mouse ATPase, H ⁺ /K ⁺ Exchanging Beta Polypeptide (ATP4b) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.066 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.