

Product Datasheet

PARP1 (phospho-T373) antibody (orb166779)

Description	Rabbit polyclonal antibody to PARP1 (phospho-T373).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	PARP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	113084
Uniprot ID	P09874
NCBI	NP_001609.2
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

