

Product Datasheet

Porcine Apolipoprotein A1 (APOA1) ELISA Kit (orb1667850)

Description	Porcine Apolipoprotein A1 (APOA1) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 28.1-1800 ng/mL. Sensitivity: 12.7 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.