

Product Datasheet

Human Adipocyte Plasma Membrane Associated Protein (APMAP) ELISA Kit (orb1667853)

Description	Human Adipocyte Plasma Membrane Associated Protein (APMAP) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.066 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.