

Product Datasheet

Rat Ankyrin Repeat Domain Protein 1 (ANKRD1) ELISA Kit (orb1667891)

Description	Rat Ankyrin Repeat Domain Protein 1 (ANKRD1) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.058 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.