

Product Datasheet

TJP2 (phospho-S978) antibody (orb166792)

Description	Rabbit polyclonal antibody to TJP2 (phospho-S978).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	TJP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	133958
Uniprot ID	Q9UDY2
NCBI	NP_001164101.1, NP_004808.2, NP_001163886.1, NP_001163887.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

