

Product Datasheet

CDKN2A (phospho-S8) antibody (orb166794)

Description	Rabbit polyclonal antibody against CDKN2A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	CDKN2A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16533
Uniprot ID	P42771
NCBI	NP_001182061.1 , NP_000068.1 , NP_478104.2
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

