

Product Datasheet

Porcine Alanine Aminotransferase (ALT) ELISA Kit (orb1667956)

Description	Porcine Alanine Aminotransferase (ALT) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 4.68-300 U/L. Sensitivity: 1.17 U/L. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.