

Product Datasheet

SQSTM1 (phospho-S403) antibody (orb166802)

Description	Rabbit polyclonal antibody to SQSTM1 (phospho-S403).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	SQSTM1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47687
Uniprot ID	Q13501
NCBI	NP_001135770.1 , NP_001135771.1 , NP_003891.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

