

Product Datasheet

General Androsterone (ADT) ELISA Kit (orb1668042)

Description General Androsterone (ADT) ELISA Kit

Reactivity All

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 1.4-90 ng/mL.
Sensitivity: 0.507 ng/mL. Assay Time:
1-4.5 h. Detection wavelength: 450
nm. Sample Volume: 50 µL. Sample
Type: serum, plasma or other
biological fluids.

Expiration Date Please enquire.