

Product Datasheet

ULK2 (phospho-S323) antibody (orb166806)

Description	Rabbit polyclonal antibody to ULK2 (phospho-S323).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	ULK2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	112694
Uniprot ID	Q8IYT8
NCBI	NP_055498.3 , NP_001136082.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

