

Product Datasheet

Porcine Angiotensin I Converting Enzyme (ACE) ELISA Kit (orb1668142)

Description	Porcine Angiotensin I Converting Enzyme (ACE) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 10.7 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.