

Product Datasheet

Human ATP Binding Cassette Transporter C10 (ABCC10) ELISA Kit (orb1668168)

Description Human ATP Binding Cassette Transporter C10 (ABCC10) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL. Sensitivity: 0.061 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.

Expiration Date Please enquire.